

# THE STANDARDIZED TRAINING CONCEPT



## Endoscopic training possibilities directly in the clinic

- The European Academy of Gynaecological Surgery (EAGS) developed the scientifically validated training program called Gynaecological Endoscopic Surgical Education and Assessment (GESEA)
- Skills in Laparoscopy and Hysteroscopy are trained to ensure patients safety and the best surgical outcome\*
- Practical skills are trained with the LASTT, SUTT and HYSTT models

**STORZ**  
KARL STORZ – ENDOSKOPE

\* Study Campo et al., Gynaecological endoscopic surgical education and assessment. A diploma programme in gynaecological endoscopic surgery

# THE STANDARDIZED TRAINING CONCEPT

## In-house Training Solutions from KARL STORZ following the GESEA Training Concept

The GESEA program makes use of three endoscopic training models that monitor the learning process and measure the skill level of an individual. LASTT provides a validated training approach for basic laparoscopic instrument handling skills like camera navigation, hand-eye coordination and bimanual coordination. SUTT focuses on laparoscopic suturing skills by offering a series of stitching and knotting exercises to train and test the quality level of needle handling, stitching, tissue approximation and knot tying. The HYSTT training model trains hysteroscopic camera navigation and instrument handling by making use of a model that simulates the female genital anatomy.



LASTT



SUTT



HYSTT

Our recommended GESEA Training Sets	
<b>CODE 1300</b>	<b>LASTT/SUTT Training Station to be used with the LASTT/SUTT Training Package, including:</b> SZABO-BERCI-SACKIER Lap. Trainer, TELE PACK X LED, TELECAM One-Chip Camera Head, Fiber Optic Light Cable, HOPKINS® Telescope 0°, HOPKINS® Telescope 30°, CLICKLINE KELLY Grasping Forceps, CLICKLINE METZENBAUM Scissors, CLICKLINE MATKOWITZ Grasping Forceps, 2x KOH Macro Needle Holders
26436	<b>LASTT Training Package</b>
26437	<b>SUTT Training Package</b>
26439	<b>SUTT1 Pads for Training</b>
26440	<b>SUTT2 Pads for Training</b>
<b>CODE 1301</b>	<b>HYSTT Training Station to be used with HYSTT Training Package, including:</b> TELE PACK X LED, TELECAM One-chip Camera Head, Fiber Optic Light Cable, HOPKINS® Telescope 30°, BETTOCCHI® Inner Sheath, BETTOCCHI® Outer Sheath, 2x Biopsy and Grasping Forceps
<b>26438</b>	<b>HYSTT Training Package</b>
<b>CODE 1311</b>	<b>Combined LASTT/SUTT and HYSTT Training Station</b>

For direct contact email us at [gyn@karlstorz.com](mailto:gyn@karlstorz.com)

It is recommended to check the suitability of the product for the intended procedure prior to use. Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.



KARL STORZ SE & Co. KG  
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
Postbox 230, 78503 Tuttlingen/Germany  
Phone: +49 7461 708-0  
Fax: +49 7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)

[www.karlstorz.com](http://www.karlstorz.com)

